

Flexibility for a fast-growing laboratory...



Uro-dipcheck 240 AUTOMATED URINE ANALYZER

Drive it the way you want !



- Throughput of 240 samples/hour
- Analyzes complete urine profile
- Continuous strip loading facility
- Large memory for 1000 results
- Automatic 'used strip' discarder
- On-board thermal printer
- Bi-directional data transfer capability



Technical Specifications¹:

Test Parameters :	pH, Nitrite, Glucose, Protein, Blood, Specific Gravity, Bilirubin, Urobilinogen, Ketone, Leucocytes
Throughput :	240 samples/hr
Principle :	Dual Wavelength reflectance photometer
Wavelengths :	550 nm, 620 nm & 720 nm
Display :	LCD Interface
Memory Capacity :	1000 samples
Printer :	On-board thermal printer
System Interface :	RS 232-C
Ambient Condition :	Temperature 10 °C - 30 °C , R. Humidity ≤ 80 %
Light source :	Super high luminosity lamp
Power :	AC 230 V \pm 10%, 50 Hzz
Power Dissipation :	45 VA
Dimensions :	394 mm (W) X 367 mm (D) X 160 mm (H) approx.
Weight :	4 kg approx.



Uro-dip 10A strips

Compatible with Uro-dipcheck 240.
 Complete Urine analysis in less than 60 seconds.
 Optimum inter-pad distance.



Ud Trol controls (Level I & II)

Ready-to-use quality control.
 Available in two levels to validate performance of urine analysis

¹Specifications are subject to change without prior notice.

ERBA diagnostics Mannheim GmbH

Mallaustrasse 69-73, D-68219 Mannheim / Germany
 Tel.: (+49) 0621 8799770 • Fax: (+49) 0621 8799688
 Email: sales@erbamannheim.com • Website: www.erbamannheim.com

